

Release notes HI6321 v1.1.1 2024-10-10

What's New

- Added “Save Bitmap” feature to allow user to save plot from Log Recall View Screen.
- Implemented new printer functionality that allows the user to select a character table to use for handling special characters.
- Enabled the use of special characters for user settings items that are not subject to restrictions
- User interface improvements:
 - Increased scroll wheels size for setup items

V1.03	V1.1.1
<div><div>Date</div><div><div>Year</div><div>2021</div><div>2022</div><div>2023</div></div><div>/</div><div><div>Month</div><div>2</div><div>3</div><div>4</div></div><div>/</div><div><div>Day</div><div>2</div><div>3</div><div>4</div></div></div>	<div><div>Year</div><div>2022</div><div>2023</div><div>2024</div><div>2025</div><div>2026</div></div> <div>/</div> <div><div>Month</div><div>7</div><div>8</div><div>9</div><div>10</div><div>11</div></div> <div>/</div> <div><div>Day</div><div>20</div><div>27</div><div>28</div><div>29</div><div>30</div></div>

- Implemented different markers styles for log types
 - Manual logs use solid circles while automatic logs use lines
- Minor updates to language translations
- Added option to hide the factory calibration warning